

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Panasil tray

Revision date: 20.04.2015

Product code: 001154

Page 1 of 4

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Panasil tray

#### Further trade names

Panasil tray fast Heavy

Panasil tray soft Heavy

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

High viscosity, addition-curing, elastomeric precision impression material based on polyvinylsiloxane. For dental use.

### 1.3. Details of the supplier of the safety data sheet

Company name:	Kettenbach GmbH & Co. KG	
Street:	Im Heerfeld 7	
Place:	D-35713 Eschenburg	
Telephone:	+ 49(0) 2774/705-0	Telefax: +49 (0)2774/705-87
e-mail:	info@kettenbach.de	
Contact person:	Regulatory Affairs Manager	Telephone: +49 (0)2774/705-0
Internet:	www.kettenbach.de	
Responsible Department:	Regulatory Affairs	

### 1.4. Emergency telephone number:

+49 (0)30 30686790 (English/German)

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Regulation (EC) No. 1272/2008

This mixture is not classified as hazardous in accordance with Regulation (EC) No. 1272/2008.

### 2.2. Label elements

#### Additional advice on labelling

There is no requirement for the product to be specially labelled according to EC directives or the corresponding national laws.

### 2.3. Other hazards

No risks worthy of mention.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Chemical characterization

Preparation of silicone polymers and fillers with platinum-catalyst (ppm-range).

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### After inhalation

not applicable

#### After contact with skin

In general the product is not skin irritative.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of troubles or persistent symptoms, consult an ophthalmologist.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Panasil tray

Revision date: 20.04.2015

Product code: 001154

Page 2 of 4

#### After ingestion

Seek medical attention if problems persist.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

##### Suitable extinguishing media

Water fog. Foam. Carbon dioxide (CO<sub>2</sub>).

#### 5.3. Advice for firefighters

Co-ordinate fire-fighting measures to the fire surroundings.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

Not required.

#### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

#### 6.3. Methods and material for containment and cleaning up

Take up mechanically.

#### 6.4. Reference to other sections

No release of hazardous substances.

### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

##### Advice on safe handling

No special precautionary measures are necessary.

##### Advice on protection against fire and explosion

No special precautionary measures are necessary.

##### Further information on handling

No special precautionary measures are necessary.

#### 7.2. Conditions for safe storage, including any incompatibilities

##### Requirements for storage rooms and vessels

No special requirements.

##### Advice on storage compatibility

Not required.

##### Further information on storage conditions

No special requirements.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

#### 8.2. Exposure controls

##### Appropriate engineering controls

not applicable

##### Protective and hygiene measures

The usual precautions when handling chemicals have to be obtained.

Avoid contact with skin, eyes and clothes.

##### Eye/face protection

Not required.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Panasil tray

Revision date: 20.04.2015

Product code: 001154

Page 3 of 4

#### Hand protection

Single-use gloves.  
Suitable material: PE (polyethylene).  
Unsuitable material: NR (Natural rubber (Caoutchouc), Natural latex).  
For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

#### Skin protection

Not required.

#### Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

#### Environmental exposure controls

Not required.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state:	Paste
Colour:	Base: orange. Catalyst: light grey.
Odour:	odourless

pH-Value:	neutral
-----------	---------

#### Changes in the physical state

Melting point:	not applicable
Initial boiling point and boiling range:	not applicable

#### Explosive properties

not Explosive.

Density (at 20 °C):	1,4 - 1,6 g/cm <sup>3</sup>
Water solubility:	insoluble

#### Test method

## SECTION 10: Stability and reactivity

### 10.4. Conditions to avoid

None by intended use.

### 10.5. Incompatible materials

None by intended use.

### 10.6. Hazardous decomposition products

None by intended use.

#### Further information

Please pay attention to the Directions for Use.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

#### Acute toxicity

Acute toxicity, oral LD50: > 2000 mg/kg Rat

#### Irritation and corrosivity

Irritant effect on the skin: Not an irritant.  
Irritant effect on the eye: Not an irritant.

#### Sensitising effects

no danger of sensitization.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Panasil tray

Revision date: 20.04.2015

Product code: 001154

Page 4 of 4

## SECTION 12: Ecological information

### 12.1. Toxicity

product has not been tested.

#### Further information

Do not allow to enter into surface water or drains.

1 (slightly hazardous to water (WGK 1))

Doesn't get into the sewage water as long as the process is carried out according to regulations.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Advice on disposal

Remove according to the regulations.

Additional information:

#### Waste disposal number of waste from residues/unused products

070217 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of plastics, synthetic rubber and man-made fibres; waste containing silicones other than those mentioned in 07 02 16

#### Waste disposal number of used product

200399 MUNICIPAL WASTES (HOUSEHOLD WASTE AND SIMILAR COMMERCIAL, INDUSTRIAL AND INSTITUTIONAL WASTES) INCLUDING SEPARATELY COLLECTED FRACTIONS; other municipal wastes; municipal wastes not otherwise specified

#### Contaminated packaging

Completely emptied packings can be re-cycled.

## SECTION 14: Transport information

### Other applicable information

No dangerous good in sense of these transport regulations.

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### National regulatory information

Water contaminating class (D): 1 - slightly water contaminating

## SECTION 16: Other information

### Further Information

The information is based on present level of our knowledge. It does not, however, give assurances of product properties and establishes no contract legal rights.

The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*